

FOR IMMEDIATE RELEASE

CONTACT: Heather Ulibarri
Product Manager
Biotrol International
1-800-822-8550

August 10, 2000

Biotrol International Introduces **Biocide_{G30}TM** **30-Day Re-Usable Sterilizing and Disinfecting Solution**

Ready-To-Use - Requires No Activation!

LOUISVILLE, CO — Biotrol International, makers of Birex_{se}[®] Surface Disinfectant, and other infection control products for the dental industry, announces the introduction of Biocide_{G30}TM 30-day Re-usable Sterilizing and Disinfecting solution.

A ready-to-use, 2.65% glutaraldehyde solution, Biocide_{G30}TM requires no activator and is effective for 30 days. Biocide_{G30}TM is easy, quick and clean to use, with no concerns about mixing or use life.

Biocide_{G30}TM Solution, a slightly acidic (pH 6.3), aqueous glutaraldehyde, contains a non-ionic surfactant, a silicone antifoam agent, a corrosion inhibitor, and a lemon scent. It is packaged in convenient gallon bottles.

When used full strength, Biocide_{G30}TM provides high-level disinfection of dental instruments in 45 minutes, and complete sterilization in 10 hours. Biocide_{G30}TM Solution is proven effective against spores, vegetative organisms, fungi, and enveloped and non-enveloped viruses. In addition, Biocide_{G30}TM Solution has passed extensive materials compatibility tests and is recommended for sterilizing dental devices containing metals, plastics and elastomers.

Biocide_{G30}TM is available through authorized Biotrol dealers and is a part of the Biotrol Infection Control Instrument Care Line. For more information about Biotrol International and its products, please call 1-800-822-8550, 7:00 a.m.-5:00 p.m. Mountain Standard Time.